<Znak: 92641385-FA>

29 lipca 2021 r.

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Z poważaniem,

Marie Pierre Barlanga  
Quality Department  
Boston Scientific International S.A.

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- Insert the white non-vented spike into the saline bag (or bottle). Hang the saline bag on the bag hook.
  - Insert the white vented spike into the empty 20mL vial.
  - Remove the RED RUBBER cap from the needle injector assembly. Place the needle injector assembly on a sterile surface.
  - Slowly fill and discharge the syringe to remove air from the Administration Set tubing and syringe. Continue priming vigorously with full pressure until there are no bubbles in the lines and there are continuous streams of saline flowing out of both needle holes in the needle injector assembly.
  - Fill the syringe when priming is complete.
- 3. Dose Vial Preparation**
- Lift the TheraSphere™ dose vial in its lead pot and tilt the lead pot back and forth to 90 degrees to wet any microspheres on the vial septum. Tap the bottom of the lead pot firmly on a hard surface. Place the lead pot into the pot holder in the acrylic box base.
  - Remove the lead pot lid and place it upside down on a non-sterile surface.
  - Use a hemostat to remove the purple seal from the top of the dose vial acrylic shield. Discard the seal in the Nalgene waste container.
  - Use a sterile adhesive strip to remove the dose vial acrylic shield plug. Discard the plug and sterile adhesive strip in the Nalgene waste container.
  - Use an alcohol swab and a hemostat to swab the dose vial septum. Discard the swab in the Nalgene waste container.
  - Record the dosimeter initial reading for the dose vial (mR/h).
  - Measure and record the initial radiation field for the patient, using an ionization survey meter.
- 4. Final Assembly**
- Close the white pinch clamp on the outlet tubing between labels 'D' and 'E'.
  - Place the empty 20 mL vial in the holder on the acrylic box and push the pressure relief valve tube into gripper clip 'A'.
  - Insert the needle injector assembly into the acrylic dose vial shield. Press on the GREEN cap to lock it in place. You will hear or feel a click or snap.
  - Place the inlet tubing through slot 'B' in the acrylic box. Place the outlet tubing through slot 'D' in the acrylic box. Loop the tubing around the side and place the fitting into the holder at 'C'.
  - Clamp the priming line at label 'C' with the blue pinch clamp. For sets with no blue pinch clamp, clamp the priming line with hemostats (or equivalent).
  - Push the YELLOW tabs on the needle injector assembly all the way down, locking the needles into the dose vial. You will hear or feel a click or snap at the bottom of travel.
  - Ensure that the side shield is installed on the acrylic box. Place the top shield on the acrylic box with the sloped shield towards slot 'D'. Ensure that the tubing is not pinched or kinked.
  - Move the cart close to the patient. Lower the bed to lowest position.
  - Place a sterile towel under the extension arm holder 'E', and under holder 'C'.
  - Place a sterile towel across the gap between the acrylic box and the patient.
  - The Interventional Radiologist (IR) will flush the infusion catheter to ensure flow. Replace the infusion catheter if it is damaged or does not have satisfactory flow. Do not use a catheter extension or extra fittings. Replace the catheter if it is too short.
  - Disconnect the outlet tubing labeled 'E' from the priming tubing at holder 'C'. Firmly connect the outlet tubing 'E' to the catheter.
  - Place the catheter connection into the slotted holder 'E' at the end of the extended arm. Outlet tubing 'E' must be above the holder, with the infusion catheter hanging vertically below.
  - The IR will verify the infusion catheter position.
  - Release the white pinch clamp from the outlet tubing. Dents in tubing may be reduced by rolling outlet tubing with fingers.
- 5. TheraSphere™ Administration**
- ATTENTION:** Beta radiation fields can be very high during microsphere transfer. Stand behind beta shielding or maintain distance.
- Record the starting time of the administration.
  - Infuse TheraSphere™ Y-90 glass microspheres using steady pressure on the syringe plunger. Infuse continuously until the syringe is empty ( $\geq 20$  cc per minute).
- NOTE:** If the infusion pressure is over 30 psi, excess fluid will drip into the vented 20 mL vial. If this oc-

## Troubleshooting

Problem	Action
1. Difficulty priming the Administration Set.	<p>Verify that the tubing in the Administration Set is not pinched or kinked. Verify that the pinch clamp is not closed.</p> <p>The first priming flush should be performed very slowly to prevent small bubbles from forming in tubing and fittings. Subsequent priming flushes should be vigorous with full pressure.</p> <p>If saline leakage is observed, ensure connections are tight.</p> <p>If the issue cannot be identified and corrected, replace the Administration Set with a new one. Notify the manufacturer of the problem.</p>
2. Leakage that may contain microspheres	<p><b>Attention:</b> Any leakage from the dose vial, injector assembly, tubing 'D' through 'E', or the catheter connection at 'E' is likely to contain microspheres.</p> <p>Assess the extent of the leak. Ensure that the needle injector is properly inserted into the dose vial. If warranted, abort the infusion, disassemble the Administration Set and commence decontamination procedures. During decontamination, investigate the cause of the leak.</p>
3. Leakage of saline during infusion	<p>Leakage observed from the syringe, the saline bag/bottle, or tubing lines 'A', 'B' and 'C' will only contain saline. If saline leakage is observed during TheraSphere™ Administration, maintain steady pressure on the syringe. <b>Do not stop the flush.</b> At the end of the flush, address the saline leakage. Ensure that priming tube 'C' is clamped. Ensure connection to the syringe is tight. Adjust the saline bag or bottle connection.</p>
4. Blood begins to flow back to the TheraSphere™ dose vial, when the catheter is connected and the syringe is not being pushed.	<p>This indicates that one of the fittings or the TheraSphere™ dose vial septum is compromised. The procedure should be aborted if the issue cannot be identified and corrected. If issue has been identified and corrected, continue with administrations and observe the system for possible leaks (see Problem 2).</p>
5. Excessive fluid flow resistance is experienced during infusion  or  Difficulty achieving the desired dosimeter reading.	<p>Verify that the white pinch clamp is open. Verify that the tubing between the syringe and dose vial are not pinched or kinked. Verify that the tubing between the dose vial and catheter are not pinched or kinked. Verify that the yellow tabs are pushed all the way down.</p> <p>Apply sufficient pressure on the syringe to cause fluid to flow into the pressure relief vial.</p> <p>Apply and release pressure on the syringe several times rapidly. This may clear a collection of microspheres at the tip of the outlet needle.</p> <p><u>Close the white pinch clamp</u> before performing any actions with the catheter. Verify that there is no blood coagulation or damage in the catheter.</p> <p><b>Attention:</b> There may be microspheres in the outlet line and catheter. Use standard radiation safety methods to assess the components before handling. Use remote handling tools as appropriate.</p>

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