**Philips Healthcare**

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-1/2-

FSN86600010

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## NEGATIVE PRESSURE LIMITING (NPL) VALVE IN SIMV & VCV MODES

This information applies to Philips anesthesia machines: Siesta i Breasy, Siesta i Whispa, Siesta i TS, Dameca MRI508 and IntelliSave AX700. The content below will be added to the next revision of the IFU:

### SECTION 1.6 BUILT-IN SAFETY FEATURES

The NPL valve opens if spontaneous breathing of the patient generates a negative pressure of -5 to -7.5 cmH<sub>2</sub>O (hPa). This opening pressure of the NPL valve is defined at an inspiratory flow rate of 3.5-4.5 L/min. Negative airway pressure is only possible if the patient takes a large spontaneous breath around the same time as a mandatory breath.

### SECTION 2.1.2 PATIENT SAFETY

If the patient is connected to the anesthesia machine in a volume controlled ventilation mode, and attempts a large spontaneous breath around the same time as a mandatory breath, it is possible for the patient to create a high negative airway pressure. If the negative airway pressure is below -7.5 cmH<sub>2</sub>O (hPa), the NPL valve of the anesthesia machine will open, but due to the resistance of the flow and the time taken to reach the patient, the airway pressure may decrease further. The negative pressure however is not sustained. If the patient is breathing spontaneously, and if the machine is installed with support ventilation modes, VSV or PSV could be considered.

### SECTION 4.3-A VENTILATION MODES (SECTION 4.3.1 IN SIESTA I BREASY)

#### NOTE

If the patient is connected to the anesthesia machine in a volume controlled ventilation mode, and attempts a large spontaneous breath around the same time as a mandatory breath, it is possible for the patient to create a high negative airway pressure. If the negative airway pressure is below -7.5 cmH<sub>2</sub>O (hPa), the NPL valve of the anesthesia machine will open, but due to the resistance of the flow and the time taken to reach the patient, the airway pressure may decrease further. The negative pressure however is not sustained. If the patient is breathing spontaneously, and if the machine is installed with support ventilation modes, VSV or PSV could be considered.

### SECTION 9.4 VENTILATOR - PRESSURE RANGE

NPL valve opening pressure: -5 to -7.5 cmH<sub>2</sub>O (hPa) at a flow rate of 3.5 to 4.5 L/min.

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© 2014 Philips Philips reserves the right to modify technical specification and design without prior notice
Philips Anesthesia Care A/S, Islevdalvej 211, 2610 Roedovre, Denmark
Dameca A/S, Islevdalvej 211, 2610 Roedovre, Denmark
Philips Anesthesia Care A/S