20.10.2016

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CE 0499

Product description:

ROTAIO is a disc replacement implant for use in the cervical spine. In addition to restoring height of the vertebral disc, its primary function is to preserve physiological mobility in the affected segment.

The prosthesis consists of a superior and an inferior endplate on which the sliding elements are anchored and secured by means of a fixation pin. The implant design provides an optimum fit to the anatomy of the intervertebral space. For optimum primary stabilisation, the endplates have a toothed surface, and special surface treatment of the endplates assists bony integration. To provide maximum coverage of the vertebral body endplates, a total of 16 different disc prostheses are available with various footprints and heights.

Instruments specially developed by SIGNUS, that ensure safe application, are available for application of the implant system.

Indications:

ROTAIO is a disc replacement implant for use in the cervical spine (C3 – C7). In addition to restoring height of the vertebral disc, this implant, in contrast to a cage, serves to preserve physiological mobility. The implant is used for the following indications following cervical discectomy via anterior access:

- Discopathy
- Disc herniation
- Foraminal and spinal canal stenosis

Contraindications:

- Instability
- Severe facet joint and spinal disc degeneration
- Immobility of the affected segment
- History of surgery to the affected segment
- Deformity
- Ossification of the posterior longitudinal ligament
- Traumatic lesions of the cervical spine
- Local tumour
- Osteoporosis, osteopenia
- Acute or chronic systemic, spinal or localized infections
- Systemic or metabolic diseases
- Allergy or intolerance to implant material (e.g. nickel)
- Surgical conditions which rule out any potential benefit from spinal surgery (such as severe damage to bone structures at the implantation site, badly distorted anatomy due to anomalies)
- Medical conditions that could prevent successful implantation (e.g., obesity, mental illness, pregnancy, paediatric cases, patients in poor general health, lack of patient compliance)
- Cases not mentioned under Indications

Material:

The disc prosthesis is made of the following materials:

- Superior and inferior endplates: Titanium alloy as per ISO 5832-3
- Sliding components: Cobalt-chrome alloy as per ISO 5832-12

The material is an established material for use as an implant. It is biocompatible, corrosion-resistant, non-toxic in the biological environment and allows for interference-free X-ray imaging. It is MRI-compatible with reservations. Safety information must be obtained from SIGNUS before an MR examination is preformed.

Sterility:

All implants are supplied in fully preassembled condition in double sterile packaging, and are gamma-sterilised in accordance with DIN EN ISO 11377 et seq.

The instruments GB0006, GB14 and GB58, which are supplied in sterile condition as part of the set, are intended for single use only and are not reusable. They can be reprocessed (e.g. after accidentally opening the packaging) with one of the commonly used sterilisation methods, provided that this is compatible with the hospital guidelines. It is essential that all packaging components are removed completely prior to resterilisation.

ROTAIO is intended for single use only and is not reusable. The implant must not be reesterilised.
Reprocessing and/ or reuse can result in infection and/ or loss of function, which can in extreme cases lead to the death of the patient.

Precautions:

- Store sterile implants and sterile instruments in their original packaging.
- Do not remove implants from their sterile, protective packing until the moment they are required for use.
- Before use, check the expiry date and the integrity of the sterile packaging.
- Prior to implantation, implants and instruments must also be checked for integrity and the given size compared against the size determined using the trials.

Reprocessing of instruments:

- Completely remove all components of the packaging prior to reprocessing
- Observe the validated reprocessing procedure in the instructions included with the tray
- The instrument tray must undergo a validated cleaning process before being returned to SIGNUS. This must be documented on the delivery note provided and enclosed with the return shipment.

Labelling:

Explanation of the symbols that may be used on SIGNUS product labels:

- CE 0499** CE mark
- Do not reuse
- Catalogue number
- Use by
- Do not sterilise
- Manufacturer
- Sterilised by irradiation
- Lot number
- Consult instructions for use
- Do not use if packaging is damaged

Application:

- The attending physician, who must be experienced in the implementation of surgical interventions and instructed in the use of the instruments, is responsible for determining the indication, selecting the implant and performing the implantation.
- All information about the surgical technique, the range of implants, the instruments and their use is described in detail in the SIGNUS product information. This information must be available on-site and must be known to the surgical team.
- Before performing the surgical intervention, it must be ensured that all required implants and instruments are available on-site and that they are in good working order.
- If there are any pre-operative uncertainties relating to the implant system, information must be obtained from SIGNUS.
- Before the surgical intervention, the patient must be informed of all possible risks and complications that can arise in connection with the intervention itself and with use of the implant.
- The correct position of the spinal implant must be verified with the aid of imaging procedures both during the intervention and at the end.
- The implant must be checked for intactness and correct size before being implanted. Footprint: Compare the footprint information on the adapter with the information indicated on the selected trial instrument. Height: observe the dots marked on the posterior endplate of the implant (1 / 2 / 3 / 4 dot(s) = 5 / 6 / 7 / 8 mm, respectively)
- The sterile-packaged ROTAIO disc prosthesis is fitted with an adapter, which facilitates insertion into the intervertebral space. This adapter must be removed after implantation and must not remain in the patient.
- If the prosthesis is to be re-positioned, it must be grasped carefully with the extraction forceps and withdrawn. For withdrawal, the distraction in the affected segment must be increased. Excessive force (e.g. by using the slotted hammer) must not under any circumstances be exerted to avoid damaging the prosthesis. Before re-mounting the implant onto the insertion adapter, the intactness of the implant must be verified.
- If the implant has to be re-mounted onto the insertion adapter, the instructions in this regard provided in the tray or in the product information must be followed.
- The article number, designation and lot number of the implant used must be documented in the patient record. All necessary data are indicated on the labels in the original packaging and must be pasted into the patient record to ensure lot traceability.
- In the postoperative phase, special care must be taken to ensure that the patient is given all the necessary information by the treating physician according to the patient's individual requirements.
- The implant must be checked for correct size prior to the intervention.
- Particular attention must be paid to the protection of nerve roots.
- After preparation, carefully inspect the intervertebral disc cavity for bone fragments.
- Excessive ablation or even complete removal of the cortical base plates and cover plates of adjacent vertebral bodies must be avoided.
- When inserting the implant, refrain from using excessive force in order to protect the adjacent vertebral bodies.
- Do not forcibly hammer the implant into place.

Risks/adverse reactions:

These instructions for use do not list the general risks associated with surgery or the complications that can arise from spinal surgery. The following are potential risks and complications which are related to the intervertebral disc prosthesis and may necessitate repeat surgery:

- Wear or breakage of implant components
- Loss of fixation, dislocation, subsidence
- Temporary or permanent development of noise
- Sensitivity to foreign bodies, allergic reactions or other local/ systemic adverse reactions to the implant materials used
- Incorrect placement
- Infection
- Vascular lesion
- Neural lesions with reversible or permanent neurological deficits or paralysis
- Fusion of the segment that underwent surgery

These risks can potentially lead to injuries of all degrees of severity to the surrounding tissue, the nerves and blood vessels, which can in extreme cases even lead to death.

USA: Federal law restricts this device to sale by or on the order of a physician.

Product Guarantee:

SIGNUS Medizintechnik GmbH guarantees that every spinal implant is manufactured with the greatest possible care and from selected materials. All processes, including manufacture, packaging and testing, are subject to continuous quality control. SIGNUS Medizintechnik GmbH has no influence on the conditions under which the spinal implant is selected and used, the diagnosis of the patient, the surgical technique or handling of the implant once it has been dispatched from the company and therefore SIGNUS Medizintechnik GmbH cannot guarantee the outcome or the absence of complications. Please inform SIGNUS immediately of any report of (possible) malfunction, quoting the article number(s) and the lot number(s).

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