

Item code	Description	Affected Lot Numbers				
287/25015	Guedel Airway 8 CM M.3	19L1324FAX				
351/5979	Mechanical Filter MINI	19M0897FAX	20E1339FAX	20H0863FAX	20H0864FAX	20H0865FAX 20J0944FAX
		20K0404FAX				
353/19004	HME for Tracheostomized Patients	20A1048FAX	20B1331FAX	20H1151FAX		
352/5877	Adult-Pediatric Electrostatic Filter HME	20C0174FAX	20G1495FAX	20H0509FAX	20H0546FAX	20H0608FAX 20H1226FAX
		20L1363FAX	20K0943FAX	20K1417FAX	20L0513FAX	20L0678FAX 20L1359FAX
		20J0449FAX	20M0063FAX	20M0062FAX	20L1360FAX	
355/5430	Infant-Pediatric Electrostatic Filter HME	20C0869FAX	20F0673FAX	20G1377FAX	20K0220FAX	20L0340FAX
287/7843	Guedel Airway MIS. 11 CM	20D0616FAX	20M0890FAX			
291/9250	Bifurcated Connector	20D1496FAX				
350/5422	Electrostatic Filter	20F0181FAX				
350/19003	Electrostatic Filter	20G0317FAX				
350/5845	Electrostatic Filter	20G0638FAX				
444SPY1314	Closed Suction System14CH FIXED X SHILEY	20G0838FAX				
351/5410	Mechanical Filter	20G1184FAX	20M0618FAX	21A0219FAX	21A0255FAX	21A0470FAX
350/5879	Electrostatic Filter	20G1554FAX	20G1555FAX	20H0054FAX	20L1331FAX	20L1332FAX 21A0791FAX
351/5878	Mechanical Filter	20G1667FAX				
351/5856	Mechanical Filter 22M-22F	20H0428FAX				
444SP02116	Closed Suction System, Replcble Cath. MDI 16CH	20H0496FAX				
444SP01314	Closed Suction System14CH FIXED TRACHEOSTOM	20H0629FAX				
350/5882	Electrostatic Filter with Catheter Mount	20J0204FAX				
355/5427	Pediatric Neonatal Electrostatic Filter HME	20J0471FAX	21A0328FAX			
444SP02514	Closed Suction System, Replcble Cath. MDI 14CH	20J0842FAX				
287/7633	Guedel Airway MIS. 8 CM	20J0976FAX				
291/7617	Suction Valve	20J1258FAX				
351/5979TC	Mechanical Filter with Connector Small, Tethered Cap	20K0014FAX				
352/5855	Adult-Pediatric Electrostatic Filter HME with Catheter Mount	20K0522FAX				
352/5867	Adult-Pediatric Electrostatic Filter HME with Connector	20K1518FAX				
444SP01016	Closed Suction System16CH FIXED	20L0028FAX				
352/5844	Adult-Pediatric Electrostatic Filter HME	20L0407FAX				
354/5876	Adult-Pediatric Mechanical Filter HME	20L0671FAX	20M0673FAX	20M0732FAX	20M0733FAX	20M1480FAX

620/5173

Connector 'L'22M/15F-15M

20M0253FAX

351/5834

Mechanical Filter 'S'

20M0554FAX